RADIO FREQUENCY INTERFERENCE CHECK – INTOXILYZER 9000

Each instrument has an RFI antenna located behind its faceplate display. Radio frequencies can interfere with breath alcohol tests so it is important to perform this check for your agency's commonly used transmission channels. During the performance of this RFI check, the .08 calibration check solution simulator shall be used in the *Proficiency Testing* mode.

Fill out the top section of the RFI Report

Agency name
Location
Instrument serial number
Your name and date

Choose common channels that your agency transmits on. Fill in Transmitter Type codes under Radio Frequency Interference Tests. Transmitter Types are as follows:

Handheld: **HH**Base: **B**Mobile units: **M**

On the RFI report under Radio Transmission Data, enter call letters, channels, frequency (MHz), and wattage of each channel.

Use the *Proficiency Testing* mode to perform the RFI checks. Failure to use *Proficiency Testing* mode may cause the instrument to lock up. When prompted to connect a proficiency solution to the IR 9000, instead connect the .08 standard check solution.

Ensure the instrument is in "Ready Mode" to begin proficiency testing.

Check simulator temperature for proper range (33.8 to 34.2°C).

To access proficiency testing:

[Options]

Enter password [ENTER]

Hidden Menu appears

Press green 'Control Testing' button

Press green 'Proficiency Testing' button

Enter Proficiency solution number into instrument, on data sheet, and logbook

Connect simulator

Enter Agency Code: (# # A # #) [NEXT ARROW] Enter Lot # (Solution #): (YY.MM.#) [NEXT ARROW]

Enter Operator's Last Name [NEXT ARROW] Enter Operator's First Name [NEXT ARROW] Enter Operator's Middle Initial [NEXT ARROW] Enter Certification Number [NEXT ARROW]

Review/Continue: Continue will start the first Air Blank/Test

The instrument will begin with two air blanks followed by nine proficiency runs, all preceded and followed by an air blank. The instrument print out will label each proficiency run as a "Cal Check," i.e. AACACACACACACACACA, where A= air blank and C= Proficiency.

Allow the instrument to run the first three calibration checks normally, then during the fourth (or so) calibration check, transmit an agency channel and continue to transmit a different channel during each subsequent calibration check.

<u>Note:</u> If the instrument stops out and gives a printout with each RFI detection, the operator must begin the Proficiency Testing mode over for each subsequent RFI check.

Enter channel numbers of each transmitter type under Channel.

(OVER)

RADIO FREQUENCY INTERFERENCE CHECK – INTOXILYZER 9000

Under 0.08 Ethanol Results, enter the first three concentration results (in the CONTROL row) prior to testing a radio transmitter. For each radio transmitter, enter either RFI (if RFI is detected for that radio transmission) or the Ethanol concentration (e.g. .077 or .081) if RFI is not detected for a radio transmitter. Submit the RFI Report to SLD.

RADIO FREQUENCY INTERFERENCE CHECK – INTOXILYZER 9000

AGENCY:				
LOCATION:		I	NSTRUMENT S/N:	
INSPECTOR:]	DATE:	
RADIO TRANSMI	SSION DATA	4		
CALL LETTERS	CHANNEL	FREQU	UENCY (MHz)	WATTAGE
DADIO EDECLIEN	ICV INTERE	EDENCI	F TESTS	
RADIO FREQUEN TRANSMITTER TYPE			E TESTS 0.08 ETHANOL RESULTS*	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS
TRANSMITTER TYPE	CHAN	NEL	0.08 ETHANOL	REMARKS

B Base

HH Handheld (Test as close to instrument as possible)

M Mobile (Test as close to instrument as possible)

In proficiency mode, initiate test. When instrument is sampling solution, activate transmitter.