














**ELECTRONIC AIDS TO DAILY LIVING (EADL)
LOAN BANK - ITEM LIST**

| Item No. | Item | Source | Description |
|----------|---|------------------|--|
| 1 |  SensiTrac Flat Pad | Ablenet | Easy secure mounting system for flat surfaces for a variety of switches like the “Big Red” and “Jelly Bean” switches. |
| 2 |  Mini Stunt Remote Car | Enabling Devices | Switch Access Stunt Car and Platform Popular mini stunt car responds to your capability switches. Car has multiple functions. Blue tray (included) fits on a wheelchair table. |
| 3 |  Switch Adapted Friendly Firefly | Enabling Devices | A switch activated pet that will appeal to a variety of clients to encourage switch access during evaluation. Plays songs and vibrates. Also has different textures, which help develop tactile skills, each activation by the switch runs the toy 10 to 15 seconds. |
| 4 |  Joystick Switch | Enabling Devices | Allows User to Control Four different devices/functions with Joystick control. This joystick is perfect for teaching directionality as well as other spatial relationships. The joystick activates four different devices depending on the direction the switch is moved. Great switch for activating remote control cars for example. Pushing the switch to the left will activate one device, moving it backwards will activate another and so forth. |
| 5 |  Adjustable Finger Tip Switch | Enabling Devices | Switch allows user with very limited movement to activate a device with one finger movement. Switch moves with client. Sensitivity and switch modes are adjustable. Adjustable sensitivity allows the user to activate a special circuit with the touch of a finger providing a timed (0 to 120 seconds), latched or momentary output. |




**SWITCH ELECTRONIC AIDS TO DAILY LIVING (EADL)
LOAN BANK - ITEM LIST**

| Item No. | Item | Source | Description |
|----------|--|------------------|--|
| 6 |  Palm Joystick Switch | Enabling Devices | Allows control of up to 4 devices with a slight rocking action. A new ergonomic approach to joysticks, held easily in palm of your hand. |
| 7 | Vibration Switch | Adaptiva-tion | Allows people with very limited motor control to use virtually any surface as a switch by tapping. Variable sensitivity. |
| 8 |  Scatir Switch | TASH | Specialized Switch that allows control with very small movements like eye-blink, eye brow movement. Very versatile. Optical Switch with auditory feedback. Includes all mounting accessories. |
| 9 |  Pneumatic Dual Switch | TASH | A Dual Switch is two single switches in one. It has a 1/8" or 3.5mm stereo plug. A sip and puff switch - to activate one switch, you "sip" and to activate the other, you "puff". Internal saliva trap allows for easy cleaning. |
| 10 |  Penta Multiple Switch | TASH | For clients with very limited finger-tip movement. Allows control of up to 5 items. Can be used for wheelchair control. Computer access with interface. ECU access with CA – 6 cable below. 9-pin D Female Connector. For people with fine motor skills and good targeting ability. The five switches are small (.3" / 1 cm diameter) and activated by pressing down on buttons. Tactile and auditory feedback. |


**SWITCH ELECTRONIC AIDS TO DAILY LIVING (EADL)
LOAN BANK - ITEM LIST**

| Item No. | Item | Source | Description |
|----------|--|---------------|---|
| 11 |  Mini Joystick with Push Multiple Switch | TASH | <p>Chin control option. Tactile and auditory feedback. Can be used for wheelchair control. Computer access with interface. ECU access with CA – 6 cable below. 9-pin D Female Connector.</p> <p>Four switches are activated by moving the joystick in four directions; the fifth is activated by pressing down on the joystick shaft. Two tips are included: one is flexible rubber and can be cut to desired size, and the other is hard plastic. Tactile and auditory feedback.</p> |
| 12 |  10 button Remote | Smarthome.com | <p>X10 10-Button Credit Card Remote Allows wireless control of up to 10 items. Buttons are fairly small. Requires Transceiver Unit and Appliance Modules</p> |
| 13 | Receiving Modules for x-10 Systems | Smarthome.com | Receiving Modules and Transceivers for small appliance control, call systems, lights, etc. |
| 14 |  Mini Relax w/ jacks | TASH | <p>Environmental Control Unit for Infrared Devices.</p> <p>The Mini Relax is a scanning infrared transmitter that controls the TV, VCR, or any other device that operates using infrared. It will store up to six functions. Complete with large and colorful graphics, it provides visual and auditory feedback and features two scan modes.</p> |
| 15 |  | TASH | <p>Advanced Environmental Control Unit. Can control up to 4 infrared devices and 10 electrical devices with a single switch. X-10 transmission requires a transceiver and X-10 compatible modules. Various scanning modes with auditory/visual feedback. Trainable with remote controls from items you want to access.</p> |

SWITCH ELECTRONIC AIDS TO DAILY LIVING (EADL) LOAN BANK - ITEM LIST

| | | | |
|----|--|--------------------|--|
| 16 | CA-6 Cable | TASH | Allows client to use a “multi” switch with a 9-pin connector to control with 5 devices that are each operated by a single switch. Can be used for ECU devices. Computer Access with Interface. |
| 17 | CA-2 Cable | TASH | Allows use of a dual switch with 2 devices that are each operated by a single switch |
| 18 | Mounting Kit for Env. Control | TASH | Assist with mounting needs |
| 19 | Wheelchair Brackets | TASH | Assist with mounting needs |
| 20 | Accutab Med Organizer  | Walter Drake | Dispense exact pills 3 times a day, 7 days a week. A simple click drops pills into your hand, 3 times a day. Plastic 8" long, 3 1/2" high. Accutab tab fits under a cupboard with all hardware included. |
| 21 | Pillbox Timer  | Walter Drake | Built-in timer will beep to remind you. Easy-to-use, push-button countdown timer can be set up to 24 hours ahead of time, and at intervals as needed. |
| 22 | Talking Photo Album  | Attainment Company | Easy-to-use speech output device that holds standard 4" x 6" photos. The Talking Photo Album includes a total of 24 pages, each with a individual message capacity of 10 seconds. Insert photos or picture/text cards into the transparent sleeves and record corresponding messages. User squeezes Play button to activate speech. May be an option for work or activity task cues. Individual must be able to target small button to use. |

**SWITCH ELECTRONIC AIDS TO DAILY LIVING (EADL)
LOAN BANK - ITEM LIST**

| | | | |
|----|--|-------------------------------|---|
| 23 | <p>Voice Cue</p>  | <p>Attainment Company</p> | <p>Voice Cue is a discrete auditory cueing device. Record up to five messages (60 seconds total recording time) and set the clock to play back messages at pre-set times. Each message can be assigned two playback times. For example, record "Time to take your medication." and set the clock to play the message at 7 a.m. and again at 6 p.m. If the message catches you off guard and you're not sure what it said, just press the "replay" button within one minute to repeat the last message. Change messages or play time anytime you want. Requires a pen tip to record, so messages will not easily be erased or mis-recorded. Features a digital clock, volume control, and non-volatile memory (the messages are not lost even if the batteries run low).</p> |
|----|--|-------------------------------|---|

You can also borrow individual items from the Basic Switch Evaluation Kit. See separate itemized list for that Kit.

Contact Julie Mehrl, MOT, OTR/L for information on Environmental Access Loan Bank Items and Loan Bank Information.

Contact Information:

Julie Mehrl, MOT, OTR/L
 DDS- Clinical Services Bureau
 Therapy Services Unit
 Ph: (505) 975-5024
 Fax: (505) 861-0613
 E-mail: julie.mehrl@state.nm.us

Helpful Links:

- www.tashinc.com
- www.enablingdevices.com
- www.ablenetinc.com
- www.smarthomeusa.com
- www.adaptivation.com
- www.enablemart.com
- www.wdrake.com