

















Basic Switch Access Evaluation Kit Item List

Item No.	Device	Source	Description
1.	 Jelly Bean Switch	Ablenet	The Jelly Bean Switch is considered a standard in the industry. AbleNet's Jelly Bean switch features a 2 1/2-inch activation surface and is recommended for clients who can access a smaller target area.
2.	 Specs Switch	Ablenet	A small switch with three different base plates, so it can be attached in a variety of ways. The strap plate and the stretchy Velcro strap provided allow people to strap it onto a wheelchair or bedrail, or into someone's palm. The Specs switch is versatile and durable.
3.	 Flexible Switch	Adaptation	The Flexible Switch is flexible similar to the Traction Pad but is a pressure activated switch like the Pal Pad. It can be activated with the slightest touch, less than 2 ounces, whether it comes from directly above or at a shallow angle.
4.	 Tilt Switch	Enabling Devices	Attach this gravity-sensitive switch to a head, arm or leg and when the body extremity is tilted, it will activate (or deactivate) a toy or device. The switch can be adjusted for tilt sensitivity and can be used for posture training. Comes complete with an adjustable Velcro strap.
5.	 TV Remote Module	Enabling Devices	Just plug a switch or switches (not included) into the box, and it will activate the buttons of the remote. For example, if your TV remote has a switch plugged into the "channel up" jack, then every time you press the external capability switch, the channel will increase by one.
6.	 Single Switch Latch and Timer	Ablenet	Use Switch Latch and Timers whenever you need to add Modes of Control to any switch accessible toy or appliance. Use with any pre-adapted device, or use the Battery Device Adapter to create you own adapted toy or appliance. All Switch Latch and Timers include three modes of control that determine the amount of time the appliance operates once the switch is activated: <ul style="list-style-type: none"> • Timed Seconds: Runs appliance from 1 to 60 seconds with each switch activation • Timed Minutes: Runs appliance from 1 to 60 minutes with each activation • Latch: First switch activation turns appliance ON; second activation turns it OFF

Basic Switch Access Evaluation Kit Item List

Item No.	Device	Source	Description
7.	 Portable Radio – Adapted	Enabling Devices	This portable AM/FM radio is adapted for switch activation with a standard switch jack. A good choice for evaluation with clients who enjoy music.
8.	 Movement Sensor Switch	Enabling Devices	This unique adjustable sensor switch encourages even the most hypoactive child to move. Just create motion in front of the sensor and it provides momentary, latched or timed output.
9.	 The Ultimate Switch	Enabling Devices	The Ultimate Switch is a versatile and multi-faceted switch. It's designed to be used by even the most profoundly challenged or neurologically impaired individuals. It requires only the slightest touch for activation (10 grams of force), yet can easily withstand the abuse of gross body movements or severe blows. 19" Gooseneck and universal clamp included.
10.	 Gooseneck Mount	Enabling Devices	Now position a switch anywhere, equipped with a flexible, extra stiff gooseneck, you can attach it to a table or wheelchair using a one-piece all-surface clamp. Or attach the switch to the 3" circular plate with dual-lock fastening strips (Included).
11.	 Switch Adapted Pet or Car	Enabling Devices	The Switch Evaluation Kit includes one switch adapted pet or car. Most pets or vehicles include sound and movement affects. A good motivation for evaluation of switch use for many individuals.
12.	 Folding Fan	Enabling Devices	Good for individuals who require extra somatosensory stimuli. Fans provide a gentle breeze and soothing auditory stimulation. Adapted for switch activation.

Basic Switch Access Evaluation Kit Item List

13.	 FreeHand – Wireless Control	Adaptiva tion	The FreeHand is Adaptivation's most popular single-button environmental control unit. The FreeHand allows for safe control of AC appliances such as a fan, lights, etc. With no cords to contend with, the FreeHand offers less distraction and increased operator safety. Various Modes available including, momentary, latched, and timed.
14.	 Tiny Technician	Enabling Devices	Pocket Technician The pocket technician takes the guesswork out of diagnosing problems with your equipment. The technician tests batteries (N, AAA, AA, A, C, D, 9-V), devices and switches.
15.		Adaptiva tion	For use with FreeHand Wireless Transmitter above. Receives signals from the FreeHand. No wires required! Just plug the receiver into the wall and the item you want to control into the receiver. Vary versatile system for evaluation. Plug in fans, radios, lights, kitchen appliances, etc.
16.	 Attendant Call Chime	Enabling Devices	Alerts Attendants Immediately! This call chime is easy-to-use and helpful, too! The user simply presses a button when assistance is needed and the receiver sounds up to 40 feet away. The call button is also adapted for use with a capability switch. Size: 4" x 4" x 1".
17.	Battery Interrupters (2)	Ablenet	
18.	Jack Size Adaptors		
19.	Battery Tester/Tools		
20.	Instruction Manuals (17)	Per item	

Basic Switch Access Evaluation Kit contents may include additions or similar alternate items without notice.

Helpful Links:

- www.ablenetinc.com
- www.adaptivation.com
- www.enablingdevices.com
- www.enablemart.com
- www.infogrip.com

Contact: Julie Mehrl, MOTR/L Ph: (505) 975-5024 E-Mail: julie.mehrl@state.nm.us